
NOTES:

- 1. ALL INTERCONNECTING PIPE AND FITTINGS SHALL BE INSTALLED BY THE CONTRACTOR.
- 2. ALL PIPING SHALL BE INSTALLED IN ACCORDANCE WITH NFPA 99, 2005 EDITION. PIPING LAYOUT IS DIAGRAMMATIC ONLY. FINAL ARRANGEMENT WILL BE SITE SPECIFIC. CONTRACTOR WILL BE RESPONSIBLE FOR FIELD VERIFICATION AND COORDINATING ACTUAL LOCATION WITH ALL OTHER TRADES.
- 3. CONTRACTOR SHALL VERIFY MOTOR VOLTAGE, PHASE AND AMP RATINGS BEFORE STARTING ELECTRICAL INSTALLATION; AND MAKE CERTAIN THE VOLTAGE SUPPLIED BY THE HOSPITAL IS THE SAME.
- 4. THE ROOM REQUIRES DEDICATED MECHANICAL VENTILATION WITH AN ALLOWABLE TEMPERATURE VARIANCE OF 40° TO 105°F. (FOR HIGHER AMBIENTS, CONSULT FACTORY)
- 5. ELECTRICAL POWER SHALL BE SUPPLIED FROM THE EQUIPMENT SYSTEM BRANCH OF THE ESSENTIAL ELECTRICAL SYSTEM (EMERGENCY POWER).
- 6. EQUIPMENT, INSTALLATION, AND TESTING OF THE MEDICAL VACUUM SYSTEM SHALL COMPLY WITH NFPA 99, 2005 EDITION, AND ALL STATE AND LOCAL CODES OR ORDINANCES.
- 7. WIRE LAG & HIGH TEMP. ALARM REMOTE CONTACTS IN THIS CONTROL SYSTEM TO MASTER ALARM SYSTEM AS REQUIRED BY NFPA 99.
- 8. ALL ACCESSORIES, PIPE AND FITTINGS, EXCEPT DISCHARGE SILENCERS, BEYOND INLET AND DISCHARGE PORTS ARE SUPPLIED AND INSTALLED BY OTHERS.
- 9. WHEN DETERMINING THE TOTAL PIPE LENGTH, ADD ALL THE STRAIGHT LENGTHS OF PIPE TOGETHER IN ADDITION TO THE NUMBER OF ELBOWS TIME THE THE EFFECTIVE PIPE LENGTH FOR THAT PIPE SIZE. (SEE EQUIVALENT PIPE LENGTH TABLE & EXAMPLE CALCULATION IN THE O&M MANUAL)
- 10. ALWAYS CONSIDER SURROUNDING OCCUPANCIES FOR ALL ROTATING EQUIPMENT. ADDITIONAL ISOLATION/VIBRATION TRANSLATION REQUIREMENTS DUE TO SYSTEM LOCATION SHOULD BE DETERMINED AND PROVIDED BY OTHERS.



TANK DRAIN CONNECTION



NO VIBRATION PADS ARE NECESSARY AT THIS LOCATION

PIPE EXHAUST OUTSIDE, TURN DOWN AND SCREEN PER NFPA 99 -3" EXHAUST SILENCER (1 SUPPLIED WITH UNIT) -DISCHARGE 3" FLANGE

-MAIN LINE VACUUM SWITCH

FROM HOSPITAL

FOR INTAKE PIPE SIZING REFER TO STANDARD PRESSURE DROP TABLES. TOTAL PRESSURE DROP SHOULD BE LESS THAN 2" Hg AT A SYSTEM VACUUM LEVEL OF 15" Hg.

SINGLE POINT POWER WIRING ACCESS

-ALARM WIRING ACCESS

MAIN VACUUM GAUGE

SOURCE SHUTOFF VALVE

3" FLANGE

HOSPITAL CONNECTION

ISOLATION VALVE

EXPANSION PUMP

FOR FUTURE

FUTURE EXPANSION-

PUMP

ETHERNET WIRING ACCESS

EXHAUST PIPE SIZE TABLE (SEE NOTE 9) PIPE LENGTH **UNIT SIZE** 51'-150' | 151'-500' 0-50' 2.50" NPT 3.00" NPT 4.00" NPT PIPE LENGTH 351'-500' UNIT SIZE 0-75' 76'-350' TX 7.5 |3.00" NPT|4.00" NPT|5.00" NP1 PIPE LENGTH 76'-250' 251'-500' UNIT SIZE 0-75' 3.00" NPT 4.00" NPT 5.00" NPT TX 8.7

EFFECTIVE PIPE LENGTHS FOR ELBOWS PIPE SIZE (IN) | 2.50" NPT | 3.00" NPT | 4.00" NPT | 5.00" NPT EFF. PIPE 6.4 7.9 10 11.9 LENGTH (FT)

> FIELD PIPING INLET FIELD PIPING OUTLET

Date: Drawn: JDY This drawing and the information contained thereon remain the property of BeaconMedaes and may not be used for other than the purpose for which it is 4107 8516 74 **INSTALLATION DIAGRAM** 6.4-8.7HP DX-EXP-TX VSD CLAW DN VAC SPC STANDARD & O2 Sheet 1 of 1 DO NOT SCALE THIS DOCUMENT **ASSURED**

(8.7HP WITH 200 GALLON RECEIVER SHOWN)

Form F-007 Rev. 01