NOTES: INTAKE 2" FLANGE 1. ALL INTERCONNECTING PIPE AND FITTINGS SHALL BE INSTALLED BY THE CONTRACTOR. -WIRE PRESSURE TRANSDUCER AND ALL PIPING SHALL BE INSTALLED IN ACCORDANCE WITH NFPA 99, 2005 EDITION. PIPING LAYOUT IS **BACKUP SWITCH** DIAGRAMMATIC ONLY. FINAL ARRANGEMENT WILL BE SITE SPECIFIC. CONTRACTOR WILL BE TO PANEL DISCHARGE-RESPONSIBLE FOR FIELD VERIFICATION AND COORDINATING ACTUAL LOCATION WITH ALL OTHER TRADES. (SEE NOTE 10) 2" FLANGE CONTRACTOR SHALL VERIFY MOTOR VOLTAGE, PHASE AND AMP RATINGS BEFORE STARTING ELECTRICAL SINGLE POINT INSTALLATION; AND MAKE CERTAIN THE VOLTAGE SUPPLIED BY THE HOSPITAL IS THE SAME. POWER WIRING DISCHARGE TO-**ACCESS VACUUM PUMP MODULE** THE ROOM REQUIRES DEDICATED MECHANICAL VENTILATION WITH AN ALLOWABLE TEMPERATURE 2" FLANGE -ALARM WIRING VARIANCE OF 40° TO 105°F. (FOR HIGHER AMBIENTS, CONSULT FACTORY) ACCESS ELECTRICAL POWER SHALL BE SUPPLIED FROM THE EQUIPMENT SYSTEM BRANCH OF THE ESSENTIAL ELECTRICAL SYSTEM (EMERGENCY POWER). ETHERNET WIRING EQUIPMENT, INSTALLATION, AND TESTING OF THE MEDICAL VACUUM SYSTEM SHALL COMPLY WITH NFPA **ACCESS** 99, 2005 EDITION, AND ALL STATE AND LOCAL CODES OR ORDINANCES. WIRE LAG PUMP RUNNING ALARM REMOTE CONTACTS IN THIS CONTROL SYSTEM TO MASTER ALARM SYSTEM AS REQUIRED BY NFPA 99. PIPE EXHAUST OUTSIDE, TURN DOWN 8. ALL ACCESSORIES, PIPE AND FITTINGS, EXCEPT DISCHARGE SILENCERS, BEYOND INLET AND DISCHARGE AND SCREEN PER NFPA 99 PORTS ARE SUPPLIED AND INSTALLED BY OTHERS. -2" EXHAUST SILENCER WHEN DETERMINING THE TOTAL PIPE LENGTH, ADD ALL THE STRAIGHT LENGTHS OF PIPE TOGETHER (1 SUPPLIED WITH UNIT) IN ADDITION TO THE NUMBER OF ELBOWS TIME THE THE EFFECTIVE PIPE LENGTH FOR THAT PIPE SIZE. (SEE EQUIVALENT PIPE LENGTH TABLE & EXAMPLE CALCULATION IN THE O&M MANUAL) 10. PUMP MODULE CONTROL WIRING IS LOW VOLTAGE WIRING (24Vdc). LOW VOLTAGE WIRING MUST BE ROUTED SEPARATELY FROM HIGH VOLTAGE WIRING (200-480Vac). SEPARATION OF LOW VOLTAGE AND HIGH VOLTAGE WIRING MUST ALSO BE MAINTAINED INSIDE THE CONTROL PANEL. 24Vdc LOW VOLTAGE FIELD WIRING TO BE A MINIMUM OF 22 GUAGE, 75 DEG C COPPER WIRE. HOSPITAL CONNECTION-11. ALL MOTOR WIRING MUST BE RUN IN SEPARATE CONDUIT FROM LOW VOLTAGE WIRING 2" FLANGE AND ENTER THE CONTROL PANEL ON THE BOTTOM RIGHT SIDE, USING THE KNOCK OUTS PROVIDED. SOURCE SHUTOFF VALVE-MAIN VACUUM GAUGE 12. ALWAYS CONSIDER SURROUNDING OCCUPANCIES FOR ALL ROTATING EQUIPMENT. ADDITIONAL ISOLATION/VIBRATION TRANSLATION REQUIREMENTS DUE TO SYSTEM MAIN LINE VACUUM SWITCH-LOCATION SHOULD BE DETERMINED AND PROVIDED BY OTHERS. FROM HOSPITAL—► FOR INTAKE PIPE SIZING REFER TO STANDARD-NO VIBRATION PADS ARE EXHAUST PIPE SIZE TABLE PRESSURE DROP TABLES. TOTAL PRESSURE **NECESSARY AT THIS POINT** (SEE NOTE 9) DROP SHOULD BE LESS THAN 2" Hg AT A PIPE LENGTH SYSTEM VACUUM LEVEL OF 15" Hg. **UNIT SIZE** 101'-300' 301'-500' 0-100' 2.50" NPT 3.00" NPT 4.00" NPT TANK DRAIN CONNECTION DRIP LEG VALVE-24 VDC ELECTRICAL EFFECTIVE PIPE LENGTHS FOR ELBOWS FIELD PIPING INLET PIPE SIZE (IN) | 2.50" NPT | 3.00" NPT | 4.00" NPT This drawing and the information contained thereon remain the property of BeaconMedaes and may not be used for other than the purpose for which it is EFF. PIPE 6.4 7.9 10 FIELD PIPING OUTLET LENGTH (FT) 4107 8516 30 **INSTALLATION DIAGRAM 5.4HP TX** (200 GALLON TANK SHOWN) CLAW VSD VAC MOD STANDARD & O2 ASSURED DO NOT SCALE THIS DOCUMENT Form F-007 Rev. 01