## NOTES: MAIN VACUUM GAUGE--MAIN LINE VACUUM SWTICH 1. ALL INTERCONNECTING PIPE AND FITTINGS SHALL BE INSTALLED BY THE CONTRACTOR. SOURCE SHUTOFF VALVE-2. ALL PIPING SHALL BE INSTALLED IN ACCORDANCE WITH NFPA 99, 2005 EDITION. PIPING LAYOUT IS FROM HOSPITAL HOSPITAL CONNECTION-DIAGRAMMATIC ONLY. FINAL ARRANGEMENT WILL BE SITE SPECIFIC. CONTRACTOR WILL BE RESPONSIBLE 2" FLANGE FOR FIELD VERIFICATION AND COORDINATING ACTUAL LOCATION WITH ALL OTHER TRADES. FOR INTAKE PIPE SIZING REFER TO STANDARD 3. CONTRACTOR SHALL VERIFY MOTOR VOLTAGE, PHASE AND AMP RATINGS BEFORE STARTING ELECTRICAL INSTALLATION; PRESSURE DROP TABLES. TOTAL PRESSURE AND MAKE CERTAIN THE VOLTAGE SUPPLIED BY THE HOSPITAL IS THE SAME. DROP SHOULD BE LESS THAN 2" Hg AT A SYSTEM VACUUM LEVEL OF 15" Hg. 4. THE ROOM REQUIRES DEDICATED MECHANICAL VENTILATION WITH AN ALLOWABLE TEMPERATURE VARIANCE OF 40° TO 105°F. (FOR HIGHER AMBIENTS, CONSULT FACTORY) 5. ELECTRICAL POWER SHALL BE SUPPLIED FROM THE EQUIPMENT SYSTEM BRANCH OF THE ESSENTIAL ELECTRICAL SYSTEM (EMERGENCY POWER). 6. EQUIPMENT, INSTALLATION, AND TESTING OF THE MEDICAL VACUUM SYSTEM SHALL COMPLY WITH NFPA 99, 2005 -SINGLE POINT POWER EDITION, AND ALL STATE AND LOCAL CODES OR ORDINANCES. WIRING ACCESS 7. WIRE LAG PUMP RUNNING ALARM REMOTE CONTACTS IN THIS CONTROL SYSTEM TO MASTER ALARM SYSTEM AS REQUIRED BY NFPA 99. 8. ALL ACCESSORIES, PIPE AND FITTINGS, EXCEPT DISCHARGE SILENCERS, BEYOND INLET AND DISCHARGE PORTS PIPE EXHAUST ARE SUPPLIED AND INSTALLED BY OTHERS ¬ OUTSIDE, TURN DOWN ALARM WIRING ACCESS AND SCREEN PER NFPA 99 9. WHEN DETERMINING THE TOTAL PIPE LENGTH, ADD ALL THE STRAIGHT LENGTHS OF PIPE TOGETHER ETHERNET WIRING ACCESS IN ADDITION TO THE NUMBER OF ELBOWS TIME THE THE EFFECTIVE PIPE LENGTH FOR THAT PIPE SIZE. (SEE EQUIVALENT PIPE LENGTH TABLE & EXAMPLE CALCULATION IN THE O&M MANUAL) ~2" EXHAUST SILENCER (1 SUPPLIED WITH UNIT) -DISCHARGE 2" FLANGE EXHAUST PIPE SIZE TABLE (SEE NOTE 9) PIPE LENGTH 0-100' 101'-300' 301'-500' TX 5.4 HP | 2.50" NPT | 3.00" NPT | 4.00" NPT EFFECTIVE PIPE LENGTH FOR ELBOWS 2.50" NPT 3.00" NPT PIPE SIZE (IN) 4.00" NPT EFF. PIPE 6.4 7.9 10 LENGTH (FT) — — — FIELD PIPING INLET ===== FIELD PIPING OUTLET DRIP LEG VALVE NO FOUNDATION OR INERTIA PAD IS REQUIRED -NO VIBRATION PADS TANK DRAIN CONNECTION (HOUSEKEEPING PAD IS OPTIONAL) 4107 8516 19 ARE NECESSARY AT THIS POINT TX CLAW VAC SPC STANDARD & BM190346 O2 ASSURED Sheet 1 of 1

omi F-007 Rev. 01