		8	1		7		6	
	NO	TES:						
	1.	. ALL INTERCONNECTING PIPE AND FITTINGS SHALL BE INSTALLED BY THE CONTRACTOR.						
D	2.	ALL PIPING SHALL BE INSTALLED IN ACCORDANCE WITH NFPA 99. PIPING LAYOUT IS DIAGRAMMATIC ONLY. FINAL ARRANGEMENT WILL BE SITE SPECIFIC. CONTRACTOR WILL BE RESPONSIBLE FOR FIELD VERIFICATION AND COORDINATING ACTUAL LOCATION WITH ALL OTHER TRADES.						
	3.	CONTRACTOR SHALL VERIFY MOTOR VOLTAGE, PHASE AND RATINGS BEFORE STARTING ELECTRICAL INSTALLATION; A MAKE CERTAIN THE VOLTAGE SUPPLIED BY THE HOSPITAL					E SAME.	
_	4.	THE ROOM REQUIRES DEDICATED MECHANICAL VENTILATION WITH AN ALLOWABLE TEMPERATURE VARIANCE OF 40° TO 105°F. (FOR HIGHER AMBIENTS, CONSULT FACTORY)						
	5.	ELECTRICAL POWER SHALL BE SUPPLIED FROM THE EQUIPMENT SYSTEM BRANCH OF THE ESSENTIAL ELECTRICAL SYSTEM (EMERGENCY POWER).						
C	6.	EQUIPMENT, INSTALLATION, AND TESTING OF THE MEDICAL VACUUM SYSTEM SHALL COMPLY WITH NFPA 99, AND ALL STATE AND LOCAL CODES OR ORDINANCES.						
	7.	WIRE LAG PUMP RUNNING ALARM REMOTE CONTACTS IN THIS CONTROL SYSTEM TO MASTER ALARM SYSTEM AS REQUIRED BY NFPA 99.						
	8.	ALL ACCESSORIES, PIPE AND FITTINGS, EXCEPT DISCHARGE S BEYOND INLET AND DISCHARGE PORTS ARE SUPPLIED AND IN					•	
	9.	WHEN DETERMINING THE TOTAL PIPE LENGTH, ADD ALL LENGTHS OF PIPE TOGETHER IN ADDITION TO THE NUM TIMES THE EFFECTIVE PIPE LENGTH FOR THAT PIPE SIZE PIPE LENGTH TABLE & EXAMPLE CALCULATION IN THE O					ELBOWS QUIVALENT	
	10.	PUMP MODULE CONTROL WIRING IS LOW VOLTAGE WIRING (24Vdc). LOW VOLTAGE WIRING MUST BE ROUTED SEPARATELY FROM HIGH VOLTAGE WIRING (200-480Vac). SEPARATION OF LOW VOLTAGE AND HIGH VOLTAGE WIRING MUST ALSO BE MAINTAINED INSIDE THE CONTROL 24Vdc LOW VOLTAGE FIELD WIRING TO BE A MINIMUM OF 22 GAUGE, 75 DEG C COPPER WIRE.						
В	11.	ALL MOTOR WIRING MUST BE RUN IN SEPARATE CONDUIT FROM LOW VOLTAGE WIRING AND ENTER THE CONTROL PANEL ON THE BOTTOM RIGHT SIDE, USING THE KNOCK OUTS PROVIDE						
		EXHAUST PIPE SIZE TABLE						
		UNIT SIZE	PIPE					
		QX 10HP	0-200' 4.00" NP1		01'-500' 00"NPT		MAIN LINE VACUU	Μ
		EFFECTIVE PIPE LENGTH FOR ELBOWS						
		PIPE SIZE (IN	ESIZE (IN) 4.00" NPT		5.00" NPT		FROM HOSPITAL	
		EFF. PIPE LENGTH (FT) 10		11.9			
А		F:	NLET		FOR INTAKE PIPE SIZING REFER TO PRESSURE DROP TABLES. TOTA DROP SHOULD BE LESS THAN SYSTEM VACUUM LEVE			
		====== F:	DUTLET	-				
		24 VDC ELECTRICAL						
L	Form F-0	007 Rev. 01 8			7		6	

