
NOTES:

- 1. ALL INTERCONNECTING PIPE AND FITTINGS SHALL BE INSTALLED BY THE CONTRACTOR.
- 2. ALL PIPING SHALL BE INSTALLED IN ACCORDANCE WITH NFPA 99. PIPING LAYOUT IS DIAGRAMMATIC ONLY. FINAL ARRANGEMENT WILL BE SITE SPECIFIC. CONTRACTOR WILL BE RESPONSIBLE FOR FIELD VERIFICATION AND COORDINATING ACTUAL LOCATION WITH ALL OTHER TRADES.
- 3. CONTRACTOR SHALL VERIFY MOTOR VOLTAGE, PHASE AND AMP RATINGS BEFORE STARTING ELECTRICAL INSTALLATION; AND MAKE CERTAIN THE VOLTAGE SUPPLIED BY THE HOSPITAL IS THE SAME.
- 4. THE ROOM REQUIRES DEDICATED MECHANICAL VENTILATION WITH AN ALLOWABLE TEMPERATURE VARIANCE OF 40° TO 105°F. (FOR HIGHER AMBIENTS, CONSULT FACTORY)
- 5. ELECTRICAL POWER SHALL BE SUPPLIED FROM THE EQUIPMENT SYSTEM BRANCH OF THE ESSENTIAL ELECTRICAL SYSTEM (EMERGENCY POWER).
- 6. EQUIPMENT, INSTALLATION, AND TESTING OF THE MEDICAL VACUUM SYSTEM SHALL COMPLY WITH NFPA 99 AND ALL STATE AND LOCAL CODES OR ORDINANCES.
- 7. WIRE LAG PUMP RUNNING ALARM REMOTE CONTACTS IN THIS CONTROL SYSTEM TO MASTER ALARM SYSTEM AS REQUIRED BY NFPA 99.
- 8. ALL ACCESSORIES, PIPE AND FITTINGS, EXCEPT DISCHARGE SILENCERS, BEYOND INLET AND DISCHARGE PORTS ARE SUPPLIED AND INSTALLED BY OTHERS.
- 9. WHEN DETERMINING THE TOTAL PIPE LENGTH, ADD ALL THE STRAIGHT LENGTHS OF PIPE TOGETHER IN ADDITION TO THE NUMBER OF ELBOWS TIMES THE EFFECTIVE PIPE LENGTH FOR THAT PIPE SIZE. (SEE EQUIVALENT PIPE LENGTH TABLE & EXAMPLE CALCULATION IN THE O&M MANUAL)
- 10. ALWAYS CONSIDER SURROUNDING OCCUPANCIES FOR ALL ROTATING EQUIPMENT. ADDITIONAL ISOLATION/VIBRATION TRANSLATION REQUIREMENTS DUE TO SYSTEM LOCATION SHOULD BE DETERMINED AND PROVIDED BY OTHERS.



NO VIBRATION PADS ARE NECESSARY AT THIS LOCATION -TANK DRAIN CONNECTION CONTRACTOR SHALL PIPE TO FLOOR DRAIN

FOR INTAKE PIPE SIZING REFER TO STANDARD-PRESSURE DROP TABLES. TOTAL PRESSURE DROP SHOULD BE LESS THAN 2" Hg AT A SYSTEM VACUUM LEVEL OF 15" Hg.

FROM HOSPITAL

EXHAUST PIPE SIZE TABLE (SEE NOTE 9) PIPE LENGTH 0-75' 76'-250' 251'-500'

TX 15 | 4.00" NPT | 5.00" NPT | 6.00" NPT

	EFFECTIVE PIPE LENGTH FOR ELBOWS				
PIPE	SIZE (IN)	4.00" NPT	5.00" NPT	6.00" NPT	
_	FF. PIPE IGTH (FT)	10	11.9	13.2	

FIELD PIPING INLET


FIELD PIPING OUTLET

ETHERNET WIRING ACCESS

ALARM WIRING ACCESS-

MAIN LINE VACUUM SWITCH -MAIN VACUUM GAUGE SOURCE SHUTOFF VALVE


-HOSPITAL CONNECTION (4" FLANGE)



SINGLE POINT POWER WIRING ACCESS

3" EXHAUST SILENCER (2 SUPPLIED WITH UNIT)

-3" DISCHARGE FLANGE (2 PLACES)



-DRIP LEG VALVE (2 PLACES)

(200 GALLON RECEIVER SHOWN)

te: Drawn: JWO This drawing and the information contained thereon remain the property of BeaconMedaes and may not be used for other than the purpose for which it is loaned without the expressed written permission 4107 8549 90 INSTALLATION DIAGRAM 15HP TX DN CLAW VSD VAC SPC STANDARD Sheet 1 of 1

Form F-007 Rev. 01