Installation, Operation and Maintenance Instructions

Alarm Zone Valve Box

Part number 4107 9027 28.01

Revision 01

April 01, 2022

Installation, Operation and Maintenance Manual

Alarm Zone Valve Box

This unit is purchased from:

Date purchased:

Model number:

Part number:

Option(s) included:

Any information, service or spare parts requests should be directed to:

BeaconMedæs
1059 Paragon Way
Rock Hill, SC 29730

Telephone: (888) 463-3427
Fax: (803) 817-5750

BeaconMedæs reserves the right to make changes and improvements to update products sold previously without notice or obligation.

Part number 4107 9027 28
Revision 01
April 01, 2022

Table of Contents

1.0 Introduction

- 1.1 Product Identification
- 1.2 User Interface Layout
 - 1.2.1 Introduction
 - 1.2.2 MEGA3 Features
 - 1.2.3 Area Alarm
 - 1.2.4 Boards
 - 1.2.5 Statement of Intended Use
- 1.3 Definition of Statements
- 1.4 Definitions
- 1.5 Configurations

2.0 Installation Procedures

- 2.1 System Installation
- 2.2 Wiring
 - 2.2.1 General Requirements
 - 2.2.2 Determining Number of Conductors
 - 2.2.3 Power Supply
- 2.3 Locally Installed Sensors

3.0 Wiring Schematics

- 3.1 NFPA Combination Wiring Diagram

4.0 Operation

- 4.1 Overview
 - 4.1.1 Main Screen
 - 4.1.2 Gas Badge
 - 4.1.3 Toolbar
 - 4.1.4 Fault Status
 - 4.1.5 Display Screens
- 4.2 Alarm Configuration
- 4.3 Setting Up an Area Alarm
 - 4.3.1 Configuring Gas Badges
 - 4.3.2 Initial Setup of an Unconfigured Area Alarm
 - 4.3.3 Miscellaneous Tab

5.0 Maintenance

6.0 Troubleshooting

Alarm Zone Valve Box

1.0 Introduction

Product Identification

Each alarm is identified by a Model number, Part number, and Lot Code.

Installation procedures vary depending on the alarm configuration.

The product ID label is located on the inside of the alarm back box (Figure 1).


Figure 1: Product Identification Labels

User Interface Layout


Figure 2: User Interface Layout

1.2.1 Introduction

The documentation found in this technical file has been compiled in order to show that the BeaconMedaes Alarm Zone Valve Box (AZVB) complies with:

- 1.NFPA 99: 2018 – all relevant requirements defined in section 5 that apply to zone valve boxes and alarms.
- 2.CSA Z7396.1: 2017 – all relevant requirements defined in section 6.
- 3.CSA C22.2 No. 205 Canadian Electrical Code, Part 2 (2nd Edition), Signal Equipment.
- 4.UL 1069, Hospital Signaling and Nurse Call Equipment, Sixth Edition (ETL)

Alarm Zone Valve Box

The documentation covers the BeaconMedaes Alarm Zone Valve Box product line listed below, ranging from a 1-valve to a 7-valve configuration. The Alarm Zone Valve Box was developed using the Atlas Copco Product Development Process. This process is designed to meet the “Design Control” requirements of ISO 13485: 2016.

ITEM NO. (ENG/SPN)	ITEM NO. (ENG/FRN)	DRAWING NO.	ITEM DESCRIPTION
4107224111 - 4116	4107224266 - 4271	4107 8571 89	ALARM ZONE VALVE BOX SINGLE...
4107224117 - 4136	4107224272 - 4291	4107 8571 90	ALARM ZONE VALVE BOX DOUBLE...
4107224137 - 4167	4107224292 - 4322	4107 8571 91	ALARM ZONE VALVE BOX TRIPLE...
4107224168 - 4190	4107224323 - 4345	4107 8571 92	ALARM ZONE VALVE BOX QUAD...
4107224191 - 4214	4107224346 - 4369	4107 8571 93	ALARM ZONE VALVE BOX FIVE...
4107224215 - 4238	4107224370 - 4393	4107 8571 94	ALARM ZONE VALVE BOX SIX...
4107224239 - 4265	4107224394 - 4420	4107 8571 95	ALARM ZONE VALVE BOX SEVEN...

*Additional part numbers applied to special configurations designed per customer request after start of production on an “as needed” basis.

BeaconMedaes Alarm Zone Valve Box Model No. Chart

AZVB - - -
 A B C

Variable	Definition	Allowable Value	Allowable Value Description
A	Size	1	Single
		2	Double
		3	Triple
		4	Quad
		5	Five
		6	Six
		7	Seven
B	Valve Sizes	A	1/2"
		B	3/4"
		C	1"
		D	1-1/4"
		E	1-1/2"
		F	2"
C	Language	ENG/SPN	English/Spanish
		ENG/FRN	English/French

Alarm Zone Valve Box

1.2.2 Alarm Zone Valve Box Features

- 7" Touch Screen LCD display.
- Red alarm MUTE INDICATOR on Screen.
- Audible alarm horn
- Heartbeat LED to indicate proper operation on internal printed circuit boards and gas sensor assemblies.

1.2.3 Area Alarm

Monitors up to 7 digital gas sensors.

Features:

- Digital pressure/vacuum read out with Low/Normal/High indicators
- Rough-in box capable of holding 7 locally mounted sensors.
- Customizable Gas ID labels with location
- General fault relay that activates on any pressure/vacuum fault condition

1.2.4 Boards

Functionality is built into the Main Display board. Wire Connections are provided in the Interconnect board and transferred to the Main Display board via a single ribbon cable.

Features:

Monitors up to 7 Digital Gas Sensor Boards.

Provides 14 normally closed dry-contact relay outputs for external monitoring.

One additional normally closed dry-contact relay output is available and opens when any fault occurs.

Heartbeat indicator to indicate proper operation on the main Display Board.

B05 Power Supply Board

Features:

- Supplies 24VDC to all circuit boards in the alarm


Figure 3: B05 Power Board

B06 Digital Gas Sensor Board

Features:

- Monitors the Pressure/Vacuum from the pipeline, and provides a gas-specific digital signal for the alarm.
- Heartbeat indicator to indicate proper operation.
- Embedded in gas sensor assembly.


Figure 4: B06 Digital Gas Sensor Board

Alarm Zone Valve Box

1.2.5 Statement of Intended Use

The Alarm Zone Valve Box is designed to combine the functionality of a zone valve box and an alarm panel into one installation. The valves provide for isolation of the medical gas and vacuum pipeline in hospitals and surgery centers. The area alarm panel is installed in the same panel and provides a means to monitor the pressures of the piped gases contained in the enclosure. These Alarm Zone Valve Boxes are installed at specified locations (per NFPA99 or CSA Z7396.1) within these facilities.

The alarm within the AZVB connects to pressure transducer(s) attached to piped gasses within the enclosure and will alert at a high/low level. In the event of an emergency, a quick access window allows the valves to be shut off to cease the flow of flammable or contaminated gas, protecting the hospital staff and patients alike. The zone valve box also serves to isolate areas in the pipeline system for prescribed maintenance as needed. Each box consists of an enclosure (installed in the rough-in drywall), a variable number of valves, an area alarm panel, and a trim/window assembly. The valve boxes are installed in accessible locations specific to each area serviced, according to NFPA or CSA code.

The Alarm Zone Valve Boxes are not intended to diagnosis, prevent, monitor, treat or alleviate a disease. Nor is it compensation for injury or handicap, replacement or modification of anatomy, control of conception, or sustaining or supporting life.

The zone valve boxes are setup based on the specification from the hospital engineers. The size of the valve is dictated by the gas flowing through the pipeline, the pressure demand needed at the outlet, and the length of pipe from outlet to the source equipment. Valves

configurations are established by the medical gas and vacuum need in each specific area. For NFPA 99: 2018 and CSA Z7396.1-17, a separate Zone Valve Box must be installed outside of each zone it controls and, at least, one per floor, one per critical-care occupancy, and one per anesthetizing location. Otherwise it is logical to install zone valves along emergency exit routes for hospital staff with as many as necessary to break the system into discrete units of an appropriate number of patients per zone.

The area alarm panel can display up to 7 gas input channels with 3 conditions per channel (Normal, High, and Low). The Alarm has 15 relay outputs to provide the customer with a warning when an alarm triggers, 1 pair controls general fault, and each gas input channel gets its own pair to alarm at high and low conditions. Each gas within the enclosure has its own badge on the display with gas identification label, location label, and current pressure.

Definition of Statements

Statements in this manual preceded by the following words are of special significance.

-  **WARNING:** Means there is a possibility of injury or death to yourself or others.
-  **CAUTION:** Means there is a possibility of damage to unit or other property.
-  **SHOCK HAZARD:** Means there is a possibility of electric shock.
-  **ATTENTION:** Means precautions for handling electrostatic sensitive devices are to be observed.

NOTE: Indicates points of particular interest for more efficient and convenient operation.

Alarm Zone Valve Box

Definitions

Area Alarm Panel

Alarm panel that monitors medical gas and vacuum systems serving a specific area.

Dry-Contact

Electrical contact isolated or unconnected from any electrical source.

Electromagnetic Compatibility (EMC)

Verification that a product meets required standards for emissions of and immunity from electromagnetic energy in its intended environment.

General Fault Relay

Single-pole single-throw dry-contact relay output. Used to activate remote alarm or building management system. Relay will open whenever ANY audible alarm on panel is in progress. General Fault Relay will deactivate only after alarm condition is corrected and alarm panel resumes normal status.

LED

Light Emitting Diode

Local Sensors

Pressure / vacuum sensors mounted inside back box.

Configurations

- Alarm comes standard unconfigured. This is to allow the customer to program their gasses specifically.
- Gas Sensors are sold separately and will install locally in the AZVB enclosure. Up to 7 sensors
- See Introduction for Model Numbering Scheme and chart.

2.0 Installation Procedures

2.1 System Installation

Valve box assemblies should be installed at specified locations per NFPA 99, Standard for Health Care Facilities or CAN/CSA-Z7396.1, Canadian Standard for Nonflammable Medical Gas Piping Systems or AS 2896, Australian Standard for Medical Gas Systems.

All valves are supplied clean for oxygen service. Care should be taken to keep debris out of valves and pipelines during installation.

Step 1: Location

Plan sufficient space to accommodate valve box assembly. Location should be free of distractions allow clear view of gas readouts and access to valves during an emergency situations.

Step 2: Secure Box to Structural Supports

Secure box mounting brackets, located on top and bottom of box, to wall structural frame supports (supplied by others). Position valve box so front edge of box is flush to 1/4" below finished wall surface. Box should be level and installed at specified height/location on wall.


Figure 5

Step 3: Connect Valves to System Piping

Piping to be connected as shown in Figure 1. Braze tube extensions to piping system, per NFPA 99, while using heat sink techniques appropriate for protection of sealing materials in both valve and gauge port plug. Unbolting and removing valve for heat protection during installation is not recommended. During brazing, joints shall be continuously purged with oil-free dry nitrogen. Slide plastic grommets as close to the center of the Zone Valve Box before brazing to avoid heat damage. Move grommets back to original location after piping has cooled.

Caution:

Excessive heat may destroy valve seals. If damage occurs, repair kit must be installed

Step 4: Gas Service Identification

Install appropriate label on each valve in box assembly. Refer to Figure 5. Labels to be provided by others.

Alarm Zone Valve Box

Step 5: Standing Pressure Test

Perform standing pressure test on each system as specified by NFPA 99, CSA-Z7396.1 or AS 2896.

Caution:

Do not install gauge before pressure test.

Step 6: Dust Cover

After system pressure test has been completed, re-install dust cover, supplied with valve box, to keep out dirt and debris during wall finishing.

Step 7: Gauge Installation (after walls are finished)

When gauge is required, remove gauge port plug from valve tube extension. Install gauge using Teflon tape or oxygen compatible sealant on gauge threads. Be certain gauge is positioned upright for ease of observation. Refer to Figure 6.

Caution:

Ensure no pressure in pipe before removing plug.

Step 8: Install Trim Frame and Door

Secure trim frame in place with eight screws provided. (Max Torque is 4 in lbs) install window by slightly bending it until it can be inserted into tracks on aluminum edge rails. Window should fit securely be captured on all sides by frame.

For emergency access, door can be removed completely by pulling ring on door.


Figure 6

Step 9: Gas Identification and Location Labels

Install appropriate gas identification label for each valve at location provided on label holder. Mark location served by each valve on the label provided.

2.2 Wiring

2.2.1 General Requirements

1. All alarms are to be powered from life safety branch of emergency power system as required by applicable standards.
2. Protect all wiring from physical damage by raceways or conduit as required by applicable standards.
3. Wire alarm panels directly to switches or sensors as required by applicable standards.
4. Wiring runs should be made with color coded wire. Record color, signal, and source of signal for each wiring lead to aid in connection of alarm finish components.
5. Avoid installing alarm panels near radio transmitters, electrical motors, or switch gear.

Alarm Zone Valve Box

2.2.2 Determining Number of Conductors

- **Relay outputs**

Two conductors are required for each signal between relay output terminals and remote device (such as a building automation system). Refer to the wiring requirements of the connecting device when selecting the appropriate type of conductors to use (i.e., wire gauge, twisted or not, shielded or not.)

Refer to **NFPA Master/Combination Wiring Diagram (Section 3.1)** for signal boards.

2.2.3 Power Supply

A power supply assembly is provided with the alarm panel. The installer must use 18AWG copper supply wiring minimum with an insulation of not less than 1/32 inch (0.8mm) thick. Power supply shall be connected to a building installed circuit breaker. Circuit breaker shall be a maximum of 15 amps and marked as disconnecting means for the alarm panel. It is recommended that the circuit breaker be in close proximity to the alarm panel and properly selected according to local and national regulations.

Holes are provided in the top and bottom of the rough-in box for main power (Figure 8).

NO OTHER HOLES SHOULD BE PUNCHED OR USED.

Holes or knockouts are provided on the top panel and bottom for entrance of low voltage field wiring (Figure 7).


Figure 7: Wire Routing

NOTE:

All wiring shall be protected from physical damage by raceways or conduit in accordance with NFPA 70, National Electric Code.

NOTE:

AC Power Wires must run thru the conduit to the metal case of prepunch hole, Connect AC power wire to AC power screw terminals, GND terminal is connected to Building Earth Directly

WARNING:

Ground wire mounting location should not be serviced while the unit is powered on

Alarm panels require 100 to 250 VAC 0/60 Hz 250mA power.

For NFPA 99 compliance, alarm panel must be connected to life safety branch of the emergency electrical system.

Alarm Zone Valve Box

1. Remove four screws from plastic power supply shield.
2. Remove plastic shield from power supply.
3. Connect incoming line and neutral wires to the terminal block. Ground wire connects directly to the rough-in box (chassis) at ground stud shown in fig. 9.
4. DO NOT remove or alter the Green factory installed ground wire from the terminal block to chassis.
5. Reinstall plastic power supply shield while making sure all high voltage wires are contained within plastic shield.
6. Secure plastic shield with four screws.


Figure 8: Supply Power Wiring

WARNING:

RISK OF ELECTRIC SHOCK

Disconnect power at the circuit breaker before removing power supply shield.

Each Alarm Zone Valve Box contains individual normally closed relay contacts for Low and High Pressure conditions for each gas, and one normally closed relay for a General Fault Condition.

Identify each pair of field installed relay output signals wires inside the alarm rough-in box. Route each pair of relay output wires as shown in reference wiring diagram image figure here to the appropriate terminals on the interconnect board. Connect each set of signal wires to the terminal blocks, noting the correct polarity(+/-). Refer to Wiring Diagram Information Here

Figure 9: Interconnect Board

2.3 Locally Installed Sensors

1. Remove sensor module from shipping carton.
2. Connect sensor with DISS fitting to the appropriate DISS check at the top of the rough-in box. Verify Gas ID labels match between the sensor and piping to ensure no cross connections occur. Repeat this process for all sensors within alarm panel.


Figure 10: Local Sensor Mounting

Alarm Zone Valve Box

NOTE:

If gas sensor wires are landed on the terminal blocks in the same order as defined by the model number of the alarm panel, then initial setup will be easier and gas sensor channels will not need to be reconfigured.

If gas sensor wires are landed on the terminal blocks from top to bottom, then initial setup will be easier.

WARNING:

Output voltage range of the sensors is 0 to 24V Max DC.

Wrap threaded end of 90° elbow with teflon tape and insert into port. Tighten 90° elbow until opening faces away from valve.

Verify the gas label on the valve piping is the same as the gas identification found on the check valve. Install Check Valve into 90° elbow, then attach sensor with DISS Connection to Check Valve.


Figure 11: Local Sensor Mounting

3.0 Wiring Schematics

3.1 NFPA Combination Wiring Diagram


Figure 12: Wiring Diagram for NFPA MEGA3 Combination Alarm

4.0 Operation

4.1 Overview

The is a touchscreen gas pressure monitoring alarm created for ease of use.


Figure 13: Main Screen

4.1.1 Main Screen

The first screen that users see on a configured unit is the MAIN screen. The MAIN screen is the main view for the alarm and shows the pressure status of all of the gases (area) the alarm is monitoring.

4.1.2 Gas Badge

Every gas monitored by the alarm will have a gas badge which will show the following information:

- Gas being monitored
- Location where gas is being used
- Gas pressure value
- Gas pressure status (Low, Normal, or High)

Figure 14: Gas Badge

4.1.3 Toolbar

Running from left to right across the bottom of the MAIN screen is the toolbar. The Normal toolbar is available when the alarm is operating in normal status and shows:


Figure 15: Normal Toolbar

To switch to the Icon toolbar, touch the HOME button.


Figure 16: Icon Toolbar

Alarm Zone Valve Box

Icons on the toolbar:

1. Home button. This icon button toggles to the MAIN screen.
2. Configuration button. This icon button toggles to the CONFIGURATION screen.
3. Information button. This icon button toggles to the INFORMATION screen.

The button that corresponds to the screen currently in use will be blue.

4.1.4 Fault Status

If a gas sensor is measuring above or below the alarm level, the high or low pressure indicator on the gas badge will flash for the gas pressure out of range.

If a gas sensor becomes disconnected or shorted, both high and low indicators for that gas will flash and an error code and message will be displayed.

The Icon toolbar will appear, and the Alarm button will be flashing red.

Figure 17: Fault Status

NOTE:

The Normal toolbar will automatically reappear when all faults are resolved.

4.1.5 Display Screens

Configuration Screen

The CONFIGURATION screen is used to configure the alarm. The configuration screen is only accessible with a password. It is divided into three parts.

The top left section is the Display section. It contains a small representative view of the MAIN screen and buttons for configuring the MAIN screen.

The top right section is the Settings section.

The bottom section is the Icon toolbar.


Figure 18: Configuration Screen

Information Screen

The INFORMATION screen provides pertinent information about the alarm – including high and low alarm values and firmware details – in one central place.

Alarm Zone Valve Box

Figure 19: Information Screen

Version Tab

Touch the VERSION tab to view details about the:

- Series, code, and style of the alarm
- Firmware version for the main display board and additional boards


Figure 20: Version Tab


NOTE:

The version screens automatically scroll. Touch the blinking up or down arrow in the lower right to stop the automatic scroll and start the manual scroll.

Touch the X button in the upper right to close and return to the INFORMATION screen.

Alarm Points Tab

Touch the ALARM POINTS tab for a quick view of the alarm points—the low and high set points—for each gas badge (area alarm) or source signal configuration (master alarm).

Area Alarm

Figure 21: Alarm Points Tab

This screen provides a complete view of the configuration and is a fast way to check items such as:

- Location where the pressure sensors are physically wired.
- Gases Monitored
- Areas Monitored
- Low and high pressure set points for area signals

Touch the X button in the upper right to close and return to the INFORMATION screen.

Test Alarm Tab

Touch the TEST ALARM tab to run a quick test of the hardware.

The audible alarm will sound for one-to-two seconds to test that the horn is working properly.

Alarm Zone Valve Box

Warnings Tab

Touch the WARNINGS tab to view any possible active warnings that would include or when a sensor input is wired, but not assigned/programmed to an alarm sign.


Figure 22: Warnings Tab

4.2 Alarm Configuration

To set up or change the configuration of an alarm, start on the CONFIGURATION screen. Users are always prompted to enter the password when accessing this screen.

Figure 23: Password entry

Type the default password (121212), and touch the Enter button. You will be prompted to change the password if not done so previously (Figure 22).

The CONFIGURATION screen displays.

Normal Status Message/Panel ID Description

The alarm allows for customization of certain labels and messages.

Normal Status Message

To change the message that is communicated when the alarm is operating in normal status:

1. Touch the [A] button that is located to the left of the configuration screen (Figure 22).
2. Touch the Erase button on the keyboard to delete the existing message (Figure 24).


Figure 24: Configuration Screen


Figure 25: Normal Status Message screen

3. Enter the new word or message.
4. Touch the Save button.

Figure 26: Panel ID Description screen

Panel ID Description

To change the description of the location of the alarm panel:

1. Touch the [B] button that is located to the right of the configuration screen (Figure 23).
2. Touch the Erase button on the keyboard to delete the existing description (Figure 25).
3. Enter the new description.
4. Touch the Save button.


Figure 27: Area Alarm User Interface

4.3 Setting Up an Area Alarm

4.3.1 Configuring Gas Badges

To configure the first badge to custom settings, touch the [1] button on the configuration screen (Figure 23).

Alarm Zone Valve Box

Eight tabs are used for configuration (Figure 27). Tabs are listed in alphabetical order.

A = [GAS INPUT]	To change which wired input should be shown on the badge
B = GAS DESC	To change color codes based on standard and country
C = [GAS NAME]	To label the gas
D = AREA MONITORED	To edit the location label for a specific gas
E = [UNITS OF MEASURE]	To change the displayed unit of measure for the pressure/vacuum
F = LOW ALARM SETTINGS	To adjust the value that should indicate that the pressure is low
G = HIGH ALARM SETTINGS	To adjust the value that should indicate that the pressure is high
H = ZERO OFFSET	To make small adjustments to the displayed pressure based on user preference


Figure 28: Configuring Gas Badges

B = Gas Desc Tab


Figure 30: Gas Desc Tab

A = [Gas Input] Tab


Figure 29: Gas Input Tab

To change which wired sensor should be displayed on that badge:

1. Touch the tab for the appropriate wired gas sensor.
2. Touch the Save button to close the screen.

To change color codes based on your standard and country:

1. Scroll through the codes using the right and left arrows at the bottom of the screen.
2. Touch the tab for the appropriate code.
3. Touch the Save button when you see a check beside your selection.

Note: If the alarm default is correct, touch the X button in the upper right to close the screen.

Alarm Zone Valve Box

C = [Gas Name] Tab

Figure 31: Gas Name Tab

Use the keyboard to change the gas name. Touch the [**<ÁÑ**] to toggle to symbols and extended Latin Characters. Touch the [**1AZ**] key to toggle back to letters and numbers.

D = Area Monitored Tab

Figure 32: Area Monitored Tab

- To add locations to be monitored:
1. Touch an unused location tab.
 2. Using the keyboard, enter the name of the location that you would like to monitor (e.g., ROOM 1). Enter multiple locations as needed (e.g., Room 1, Room 2, Room 5).
 3. A check appears when the location is selected, the save button has to be pressed if the selection is new.

To edit the location information, touch the Pencil button.

E = [Units of Measure] Tab

Figure 33: Units of Measure Tab

To change the unit of measure:

1. Touch the tab for the appropriate unit of measure.
2. Touch the Save button to close the screen.

F = Low Alarm Settings Tab


Figure 34: Low Alarm Settings Tab

You can remove the low alarm set point or customize the low alarm set point.

To remove the low alarm set point:

1. Touch the NONE tab.
2. Touch the check on the CONFIRM SETTINGS CHANGE screen.

To customize the low alarm set point:

1. Touch the ADJUST tab.
2. Enter the new set point.
3. Touch the Save button.

Alarm Zone Valve Box

G = High Alarm Settings Tab


Figure 35: High Alarm Settings Tab

You can remove the high alarm set point or customize the high alarm set point.

To remove the high alarm set point:

1. Touch the NONE tab.
2. Touch the check on the CONFIRM SETTINGS CHANGE screen.

To customize the high alarm set point:

1. Touch the ADJUST tab.
2. Enter the new set point.
3. Touch the Save button.

H = Zero Offset Tab


Figure 36: Zero Offset Tab

In some cases, facility dial gauges may show slightly different readings than the alarm (e.g., 46 PSI versus 44 PSI). To match the alarm exactly to a facility dial gauge:

1. Select the offset amount.
2. Touch the Save button to close the screen.

Configuring Additional Gas Badges

Configure additional gas badges by touching the number on the main configuration screen that corresponds with the gas badge and repeating the steps listed above.


Alarm Zone Valve Box

4.3.2 Initial Setup for an Unconfigured Area Alarm

The SETUP INFORMATION screen is the first screen that users will see on an unconfigured alarm. (The MAIN screen is the first screen that users will see on a configured alarm.)

Password Setup - Unconfigured Alarm

- Touch the Configuration button.

Configuration Button

Figure 37: Setup Information Screen

Entering the Password

The ENTER SECURITY PASSWORD screen appears.

Figure 38: Entering the Security Password

- Enter the security password. (The default password is 121212.)
- Touch the Enter button.
- Touch Continue to advance to the CONFIGURATION screen.

Note:

If you would like to change the password, touch Change Password. Enter a 4-10 digit password, and touch the Enter button. Confirm the new password. Touch the X button in the upper right to close the SECURITY PASSWORD screen.


Figure 39: Changing the Security Password

- The CONFIGURATION screen will display as unconfigured and ready for setup.


Figure 40: Unconfigured Configuration Screen

Alarm Zone Valve Box

Primary Settings Tab

The PRIMARY SETTINGS tab allows you to select the gas inputs.

- Touch the PRIMARY SETTINGS tab


Figure 41: Alarm Type Tab

Gas Sensors Tab

The PRIMARY SETTINGS screen will now be displayed with the button to set up GAS SENSORS.


Figure 42: Gas Sensors Tab

- Touch the GAS SENSORS tab.


Figure 43: Gas Input Screen

Users may select up to 7 different gases on the GAS INPUT screen. These gases correspond to the digital gas sensors wired into the barrier blocks.

- To select the first gas, touch the D01 = NONE/BLANK tab.

Figure 44: Gas Sensors Screen

- Touch the right and left arrows at the bottom of the screen to scroll through additional gas choices, including options to create custom mixes. Select the gas that is wired into the barrier blocks for the top valve.

Alarm Zone Valve Box

Figure 45: Gas Choices

- Touch the Save button when you see a check beside your selection.
- Repeat to add additional gases.
- Touch the X button in the upper right to close the GAS SENSORS screen and return to the PRIMARY SETTINGS screen. The alarm will retain your information.
- Touch the X button in the upper right to close the PRIMARY SETTINGS screen. The alarm will retain your information.

4.3.3 Miscellaneous Tab

Users can customize certain aspects of the MEGA3 Alarm's appearance and functionality. Aspects of functionality that can be adjusted include the date and time settings.

- Touch the MISCELLANEOUS tab.


Figure 46: Miscellaneous Screen

5.0 Maintenance


Figure 47: Maintenance

Alarm Zone Valve Box

***Table 1: Gas Sensor Module Options**

Part Number	Gas Sensor Assembly
4107 4020 93	O ₂ Sensor Assembly
4107 4020 94	N ₂ O Sensor Assembly
4107 4020 95	MedAir Sensor Assembly
4107 4020 96	Vac Sensor Assembly
4107 4020 97	N ₂ Sensor Assembly
4107 4020 98	IAIR Sensor Assembly
4107 4020 99	WAGD Sensor Assembly
4107 4021 00	CO ₂ Sensor Assembly
4107 4021 01	CO ₂ -O ₂ Sensor Assembly
4107 4021 02	O ₂ -CO ₂ Sensor Assembly
4107 4021 03	HE-O ₂ Sensor Assembly
4107 4021 04	O ₂ -HE Sensor Assembly
4107 4021 05	HE Sensor Assembly
4107 4021 06	LAB-AIR Sensor Assembly
4107 4021 07	LAB-VAC Sensor Assembly
4107 4021 08	O ₂ -100 Sensor Assembly
4107 4021 09	MedAir-100 Sensor Assemb
4107 4021 10	CO ₂ -100 Sensor Assembly
4107 4021 11	AR Sensor Assembly
4107 4021 12	AGSS Sensor Assembly
4107 4021 13	O ₂ -N ₂ O Sensor Assembly
4107 4021 14	Mixed Gas

Table 3: Replacement Window Options

Part Number	Window Service Kits
4107 4021 77	AZVB1 Window Kit w/ Display
4107 4021 78	AZVB2 Window Kit w/ Display
4107 4021 79	AZVB3 Window Kit w/ Display
4107 4021 80	AZVB4 Window Kit w/ Display
4107 4021 81	AZVB5 Window Kit w/ Display
4107 4021 82	AZVB6 Window Kit w/ Display
4107 4021 83	AZVB7 Window Kit w/ Display
4107 4021 84	AZVB1 Window Kit w/o Display
4107 4021 85	AZVB2 Window Kit w/o Display
4107 4021 86	AZVB3 Window Kit w/o Display
4107 4021 87	AZVB4 Window Kit w/o Display
4107 4021 88	AZVB5 Window Kit w/o Display
4107 4021 89	AZVB6 Window Kit w/o Display
4107 4021 90	AZVB7 Window Kit w/o Display

****Table 2: I/O Boards**

Part Number	I/O Boards	Description
4107 4016 37	B05 Power Module	Power supply

Cleaning:

Use clean dry microfiber cloth or soft lint-free cloth to remove any smudges located on the touchscreen. Make sure not to apply too much pressure on the screen when cleaning it. Never use paper towels or tissue paper which contain wood fibers that will scratch any plastic surface.

Periodic Maintenance:

Inspect the alarm periodically and perform a self test (see Section 4.1.5 - Information Screen). This ensures proper operation of the alarm. Any damaged components inside the alarm should be replaced by a authorized BeaconMedæx service technician.

Window Replacement/Service:

If the Alarm Zone Valve box window display becomes damaged/non-functional but the rest of the alarm zone valve box is competent; part numbers 4107 4021 77-83 are satisfactory replacement kits.

If the Alarm Zone Valve box window outer-frame only becomes damaged but the electronic display and the rest of the alarm zone valve box is competent; part numbers 4107 4021 84-90 are satisfactory replacement kits.

6.0 Troubleshooting

Symptom	Possible Cause	Corrective Action
No LED power indicator illuminated on alarm front panel	AC power is not turned on	Check AC power source.
	Blown power supply fuse	Replace fuse.
	AC power wiring is not connected	Check AC entrance wiring at power supply terminals.
	Power supply DC wire harness is not connected properly between the power supply and the main board	Check connections on each end including orientation of cable connectors. Reconnect if necessary.
	Faulty power supply assembly	Replace power supply
No green "heartbeat" power indicator is illuminated on gas sensor	Sensor not connected to the Gas Input Board.	Verify correct sensor connection to Gas Input Board.
	Gas Input Board does not have power (no "heartbeat")	Verify that Gas Input Board has been initialized and that dip switches are set correctly. See item 2 above.
Fast green "heartbeat" power indicator is illuminated on gas sensor (2 beats per second)	Sensor is reading out of range	Verify gas pressure in pipeline. If normal, replace sensor.
Sensor active and Unassigned Warning	Sensor is wired to the sensor input board, but not assigned/programmed	Assign the wired gas sensor to a valid gas input.
F1 displayed on Area badge	Gas sensor is not connected or has shorted.	Check gas sensor connection. If connection is good, replace sensor.
F2 displayed on Area badge	Gas sensor is reading below range (fault low).	Verify gas sensor connections and pipeline pressure. If connection and pressure are good, replace sensor.
F3 displayed on Area Badge	Gas sensor is reading above range (fault high).	Verify gas sensor connections and pipeline pressure. If connection and pressure are good, replace sensor.
F4 displayed on Area Badge	Gas sensor connected is wrong type of gas.	Connect the correct type gas sensor.

Alarm Zone Valve Box

Troubleshooting (continued)

Symptom	Possible Cause	Corrective Action
Area Display pressure/ vacuum reading not matching pipeline.	Bad gauge on pipeline	Verify gauge used on pipeline
	Incorrect unit of measure	Check unit of measure, if incorrect, reprogram unit of measure.
	Gas sensor out of calibration	Adjust Zero Offset to correct reading. (Note: there is a max of 3 psi / 3 "hg offset adjustment). If sensor cannot be adjusted enough, replace sensor.
High or Low alarm, with normal pressure displaying on Area display	Alarm set points have been changed.	Verify / correct alarm set points.
Audible alarm reactivates a short period of time after Silence button has been pressed	A new alarm condition has occurred.	If alarm condition previously silenced, self corrects and then reactivates again, the alarm panel treats this as a new alarm condition.

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