

DiamondCare® Medical Gas Outlet Conversion Kit

SPECIFICATION

DiamondCare® Medical Gas Outlet Conversion Kit

The DiamondCare® Medical Gas Outlet Conversion Kits update existing Diamond®/DISS III wall and console rough-in assemblies to accept DiamondCare®, Diamond®, DISS Key-Style finish assemblies. The DiamondCare® finish assembly must be ordered separately.

Each conversion kit consists of a metal transition plate with attachment screws, and internal replacement rough-in components. Separate kits must be specified for recessed wall outlet and console outlet conversions. Kits should also be specified to convert pressure gases with secondary check or vacuum without secondary check.

Ordering Information

DiamondCare® Medical Gas Outlet Conversion Kits are ordered as individual assemblies. Each assembly is individually packaged.

Conversion Kit	Product Number
DiamondCare® Pressure Recessed Wall Outlet Conversion Kit	☐ 4107 2011 99
DiamondCare® Vacuum Recessed Wall Outlet Conversion Kit	☐ 4107 2012 00
DiamondCare® Console Wall Outlet Conversion Kit	☐ 4107 2012 01
DiamondCare® Vacuum Console Wall Outlet Conversion Kit	☐ 4107 2012 02

Figure 1: DiamondCare® Pressure Recessed Wall Outlet Conversion Kit - P/N 4107 2011 99






Figure 2: DiamondCare® Vacuum Recessed Wall Outlet Conversion Kit - P/N 4107 2012 00



**Figure 3: DiamondCare[®] Pressure Console Outlet
Conversion Kit - P/N 4107 2012 01**



**Figure 4: DiamondCare[®] Vacuum Console Outlet
Conversion Kit - P/N 4107 2012 02**

