

This drawing and the information contained thereon remain the property of BeaconMedics and may not be used for other than the purpose for which it is loaned without the expressed written permission from BeaconMedics Engineering.

Notes:







1. SYSTEM WEIGHT: 1101 LBS.
2.  DENOTES CENTER OF GRAVITY.

TABLE	
ID	DESCRIPTION
I	INLET LOCATION
D	DISCHARGE LOCATION
T	TANK DRAIN



INSPECTION CHECKLIST
UNLESS OTHERWISE SPECIFIED:

ALL ITEMS MARKED BY ONE OF THE ABOVE SYMBOLS ARE TO BE INSPECTED
EQUIPMENT TO BE USED DURING RECEIVING INSPECTION SHALL BE PER QUALITY SYSTEM PROCEDURES

Revised: EYP	Date: 12/9/2021	Drawn: jeb	Date: 4/18/2011
UNLESS OTHERWISE SPECIFIED ALL DIMENSIONS ARE IN INCHES AND TOLERANCES ARE: XXX ±.02 XXX ±.005		Scale: 4/18/2011	
Description: GENERAL ASSEMBLY DX 5HP LUBE VAC TANK MOUNTED 200G RECEIVER		Drawing Number: 4107 8530 79	
DN: HOP230756		Rev: 04	
Sheet 1 of 1		DO NOT SCALE THIS DOCUMENT	