



This drawing and the information contained thereon remain the property of BeaconMedics and may not be used for other than the purpose for which it is loaned without the expressed written permission from BeaconMedics Engineering.

Notes:

1. SYSTEM WEIGHT: 658 LBS
2.  INDICATES LIFTING POINTS
3.  DENOTES CENTER OF GRAVITY

| SYSTEM CONNECTIONS |                    |
|--------------------|--------------------|
| ID                 | DESCRIPTION        |
| I                  | INLET LOCATION     |
| D                  | DISCHARGE LOCATION |
| T                  | TANK DRAIN         |

Ø0.75 BOLTDOWNS  
(4 PLACES)

MINIMUM CLEARANCE FOR  
MAINTENANCE & VENTILATION

DISCHARGE - 1 1/4" NPT (FEMALE)  
(2 PLACES)  
(NOTE: DISCHARGE FLEX  
HOSES SHIPPED LOOSE)

HEPA INTAKE FILTER  
PER PUMP

BIOHAZARD DRAIN FLASK  
PER PUMP

HOSPITAL  
INLET CONNECTION  
1 1/2" NPT (FEMALE)



INTEGRATED  
LIFTING EYES  
SINGLE POINT  
ELECTRICAL  
CONNECTION

CONTROL CABINET

30 GALLON  
RECEIVER

RECEIVER  
BYPASS VALVE

TANK DRAIN  
CONNECTION  
1/2" NPT (FEMALE)

INSPECTION CHECKLIST  
UNLESS OTHERWISE SPECIFIED:  
   
ALL ITEMS MARKED BY ONE OF THE  
ABOVE SYMBOLS ARE TO BE INSPECTED  
EQUIPMENT TO BE USED DURING  
RECEIVING INSPECTION SHALL BE PER  
QUALITY SYSTEM PROCEDURES

|   |                 |  |                 |
|---|-----------------|--|-----------------|
| Revised: APS  | Date: 1/17/2023 | Drawn: JOL                                 | Date: 7/12/2015 |
| UNLESS OTHERWISE SPECIFIED ALL DIMENSIONS ARE IN INCHES AND TOLERANCES ARE:<br>XXX ±.02<br>XXX±.005 |                 | Scale:<br>D                                |                 |
| Description:<br>GENERAL ASSEMBLY DX 2HP LUBE<br>VTM V 30G TAE                                       |                 | Part Number:<br>4107 8553 38<br>Rev:<br>04 |                 |
| Sheet 1 of 1  |                 | DO NOT SCALE THIS DOCUMENT                 |                 |